



Contraindications & Cautions

If you have any medical or health concerns, we recommend that you consult with your doctor or healthcare provider prior to using any electrotherapy device.

Device users should always read the user instructions before use, and only use as directed or advised by a qualified medical professional. The user manual is supplied within the box packaging with each device.

Contraindications & Cautions:

The use of electrotherapy is not suitable and should NOT be used by individuals if any of the following contraindications apply;

- Anyone with a cardiac pacemaker, implantable defibrillator or any other implantable electrical device such as a cochlear hearing device.
- Anyone with a suspected or diagnosed heart condition, heart disease or cardiac disorders.
- Not suitable for babies, infants or young children.
- Should not be used on areas with active thrombosis including DVT.

The use of electrotherapy may not be suitable for individuals with any of the following conditions, and should only be used under the recommendation and strict guidance of a relevant medical professional or specialist:

- **Suspected or confirmed pregnancy.**
The electrode pads should never be placed close the abdomen or pelvis at any stage during pregnancy.
- **Suspected or diagnosed cancers of any type.**
The electrode pads should not be placed directly over areas where cancer is suspected or known to be present.
- **If you have diagnosed or suspected epilepsy or have experienced seizures.**
The electrode pads should not be placed anywhere on the head or neck, including the base of the neck.

Additional precautions should also be followed with consideration given to the treatment type, intensity level and pad placement locations.

Pad Placement Guidance:

The electrode pads can be placed safely in most areas around the body, but they should NOT be placed on any of the following areas:

- Directly over the heart.
- On the front or sides of the neck, including over the carotid arteries or jugular veins.
- On areas of the head and face that include the temples, eyes, ears, nose and inside the mouth.
- Directly on the genitals.
- Directly over broken or damaged skin.

NuroKor Devices

mibody & RX (Prescription) Mibody:

The Mibody & RX (Prescription) Mibody are fully regulated, FDA cleared and CE marked medical devices designed to deliver high quality, medically regulated products for personalised pain relief.

These products can be used from the 3rd trimester of pregnancy or for labour pain. However, the electrode pads should not be placed close to the abdomen or pelvis at any stage during pregnancy.

mitouch

The mitouch is intended for use as a personal therapy, recovery and physical wellness device. It is used to help promote recovery, support rehabilitation, manage pain and assist with performance and strength gains.

Microcurrent (MC²) programme is not recommended for use at any time during pregnancy. All other modes can be used for gentle massage benefits and low intensity pain relief from the 3rd trimester during pregnancy or for labour pain. The electrode pads should not be administered close to the abdomen or pelvis at any stage during pregnancy.

Microcurrent (MC²) is also contraindicated for anyone with suspected or diagnosed cancer of any type.